MACHEREY-NAGEL Compact photometer PF-3



- Intuitive menu guidance
 - USB connection and free software
- Add new tests anytime for free
- Lightweight and robust





MACHEREY-NAGEL www.mn-net.com

Addition to the MN family

The compact photometer PF-3 is the smallest member of the MACHEREY-NAGEL photometer family and perfectly fits our tradition of reliability, user friendliness and innovation.

The following diagram highlights the available versions and respective parameters.



The PF-3 is equipped with three LEDs emitting three different wavelengths, which are independent for each version. These versions are designed to meet the analysis requirements of specific applications. Currently, we can provide five distinct PF-3 versions for pool, drinking water, soil, COD analysis and fishing waters. Each version can analyze a particular set of parameters, which are vital to analyze in the respective area of application. All versions come in three different options including a box with the photometer, a rugged case with photometer and accessories, or a case with pre-packed reagents.



Test platforms for the PF-3

The PF-3 can handle both our convenient NANOCOLOR $^{\circledast}$ tube tests, as well as the economical VISOCOLOR $^{\circledast}$ test kits.

Combining both platforms within, the instrument provides maximum flexibility and allows customers to address their individual analysis needs.

NANOCOLOR® tube tests

Our barcoded tube tests are made of premium glass ware with an extra wide 16 mm tube opening. The accurately pre-dosed reagents cover the range of a wide variety of water and waste water parameters. The reaction methods are based on international standards (EPA, APHA, ISO, ...) and provide a large portfolio with multiple measuring ranges, that can be run on all our photometers.

- Convenient handling
- Highest accuracy and safety
- Precise and reliable results

VISOCOLOR® ECO test kits

Our *VISOCOLOR®* ECO test kits come as colorimetric or titrimetric tests and provide the basic range of water analysis tests in our portfolio. These high quality tests are packed with all necessary reagents as well as premium color charts in only one box for easy ordering. Fast dissolving reagents and the easy dosage save time and allow easy working procedures. Most colorimetric kits can be read photometrically on our compact photometers PF-12^{*Plus*} and PF-3.

- Outstanding cost-value ratio
- Time saving
- Effortless operations

VISOCOLOR® Powder Pillows

VISOCOLOR® Powder Pillows are photometric tests that combine easiest dosing of reagents with photometric precision. Each *VISOCOLOR®* Powder Pillow contains the exactly needed amount of reagents for a determination. *VISOCOLOR®* Powder Pillows can be evaluated on the compact photometers PF-12^{Plus} and PF-3.

- Photometric precision
- For MN and competitor's instruments
- Long shelf life









Customer benefits at a glance

The PF-3 is designed with the customer in mind, providing a reliable and rugged solution for water analysis right at the point of interest. While being a tough compact photometer, the PF-3 is nevertheless packed to the brim with clever technical details, making it the most advanced instrument in its class.

At the core of the PF-3 are premium high precision optics with especially selected LEDs and first-rate interference filters. Additionally, the PF-3 impresses with intuitive handling, a fully developed menu structure, as well as various power supply options and free of charge software



Small and compact

- Handy and easy to carry
- Light and mobile
- Compact dimensions

Ideal for all your analysis needs in the lab and in the field



Easy and convenient

- Intuitive operation using just 4 buttons
- Fully developed menu structure
- Bright display
- Open cuvette slot

Convenient operation without the need for time consuming training



Strong and tough

- Glass-fiber reinforced housing
- Water proof according to IP 68
- Shock-resistant optics

Safe results even under the most difficult conditions



Comfortable and fast

- Intuitive instructions with pictograms
- Color-coded reagent bottles with clear dosing instructions
- Fast-dissolving reagents save time
- Maximum safety for the user and easy disposal

Convenient handling for quick analysis

Good to know

Versatile reagent case solutions for the PF-3 allow individual combination with test kits.



Safe and flexible

- Runs on standard or rechargeable batteries
- Power supply via USB port
- Optional internal accu-pack chargeable via USB

Consistent and safe workflow



Individual and adaptive

- Various case solutions including reagents
- Rugged cases with premium foam inlays
- Flexible case content for multiple applications

Convenient operation with the perfect case for your need



Smart and clever

- Evaluation of VISOCOLOR[®] ECO tests, VISOCOLOR[®] Powder Pillows and NANOCOLOR[®] tube tests
- Results in different units
- Optional storage for just one zero-measurement
- Data transfer via USB

High measurement safety and easy data administration



Fair and economical

- Free of charge data export software
- Add new tests and parameters anytime
- Cost efficient photometer and reagent bundles
- Create customized methods for individual user

Consistently fair concept with impressive price-performance ratio

For swimming pools

Swimming pools are important areas of recreation and relaxation. To constantly provide highest water quality and a safe user environment, the water within swimming pools has to be treated, conditioned and closely monitored. For reliable and easy water analysis for swimming pools, the PF-3 is the perfect tool.

The built-in LEDs allow quality control personnel to analyze the classic swimming pool parameters chlorine, pH, cyanuric acid and total alkalinity, as well as a selection of additional parameters, thus providing fundamental and vital information on pool water quality.



Tests suitable for the PF-3 Pool

| Tests | Measurement range | Test number | No. of tests | REF |
|---------------------------------------|--|-------------|--------------|--------------------|
| VISOCOLOR [®] ECO tests | | | | |
| Alkalinity TA | 5–250 mg/L CaCO ₃ | 5-04 | 100 | 931204 |
| Bromine | 0.10–13.00 mg/L Br ₂ | 5-11 | 200 | 931211 |
| Chlorine 1, free, total and combined | 0.05–2.00 mg/L Cl ₂ | 5-35 | 150 | 931235 |
| Chlorine 2, free, total and combined* | 0.05–2.00 mg/L Cl ₂ | 5-15 | 150 | 931215 |
| free Chlorine 2* | 0.05–2.00 mg/L Cl ₂ | 5-16 | 150 | 931216 |
| Chlorine 6, free, total and combined | 0.05–6.00 mg/L Cl ₂ | 5-17 | 200 | 931217 |
| free Chlorine 6 | 0.05–6.00 mg/L Cl ₂ | 5-19 | 400 | 931219 |
| Chlorine dioxide* | 0.10–3.80 mg/L ClO ₂ | 5-21 | 150 | 931221 |
| Cyanuric acid * | 10–100 mg/L Cya | 5-23 | 100 | 931223 |
| Fluoride* | 0.1–2.0 mg/L F [−] | 5-27 | 150 | 931227 |
| Iron 1* | 0.04–2.00 mg/L Fe | 5-25 | 150 | 931225 |
| Iron 2 | 0.04–2.00 mg/L Fe | 5-26 | 150 | 931226 |
| ■ pH 6.0–8.2 | 6.10–8.40 | 5-70 | 150 | 931270 |
| Silica HR 200* | 10–200 mg/L SiO ₂ | 5-34 | 100 | 931234 |
| /ISOCOLOR [®] Powder Pillows | | | | |
| Chlorine, free | 0.03–6.00 mg/L Cl ₂ | 7-20 | 100 1000 | 936220 936220.1 |
| Chlorine, total | 0.03–6.00 mg/L Cl ₂ | 7-21 | 100 1000 | 936221 936221.1 |
| PH | pH 6.2–8.2 | 7-22 | 100 | 936222 |
| Silica LR* | 0.02–2.10 mg/L SiO ₂ | 7-24 | 100 | 936224 |
| Silica HR* | 2–210 mg/L SiO ₂ | 7-25 | 100 | 936225 |
| Sulfate* | 15–200 mg/L SO ₄ ^{2–} | 7-23 | 100 | 936223 |
| VANOCOLOR [®] tests | | | | |
| Chlorine/Ozone | 0.05–2.50 mg/L Cl_2 / 0.05–2.00 mg/L O_3 | 0-17 | 20 | 985017 |
| ■ pH 6.5–8.2 | 6.10-8.40 | 0-72 | 100 | 91872 |



Ordering information

| REF |
|--------|
| 919340 |
| 934102 |
| 934118 |
| 934119 |
| |

For drinking water

Drinking water is the most important foodstuff and cannot be replaced. Hence, it has to be wholesome and free of impurities that may adversely affect human health. Continuous water analysis is therefore the fundament for a sustainable and reliable drinking water supply.

The compact photometer PF-3 in combination with the proven *VISOCOLOR*[®] and *NANOCOLOR*[®] tests allows a quick and easy way of monitoring and checking the most important drinking water parameters (chlorine, chlorine dioxide, pH, fluoride and iron) for all water suppliers and industrial users during drinking water treatment.



Tests suitable for the PF-3 Drinking Water

| Tests | Measurement range | Test number | No. of tests | REF |
|---------------------------------------|--|-------------|--------------|--------------------|
| VISOCOLOR [®] ECO tests | | | | |
| Alkalinity TA | 5–250 mg/L CaCO ₃ | 5-04 | 100 | 931204 |
| Bromine | 0.10–13.00 mg/L Br ₂ | 5-11 | 200 | 931211 |
| Chlorine 1, free, total and combined | 0.05–2.00 mg/L Cl ₂ | 5-35 | 150 | 931235 |
| Chlorine 2, free, total and combined* | 0.05–2.00 mg/L Cl ₂ | 5-15 | 150 | 931215 |
| ■ free Chlorine 2* | 0.05–2.00 mg/L Cl ₂ | 5-16 | 150 | 931216 |
| Chlorine 6, free, total and combined | 0.05–6.00 mg/L Cl ₂ | 5-17 | 200 | 931217 |
| free Chlorine 6 | 0.05–6.00 mg/L Cl ₂ | 5-19 | 400 | 931219 |
| Chlorine dioxide* | 0.10–3.80 mg/L CIO ₂ | 5-21 | 150 | 931221 |
| Cyanuric acid * | 10–100 mg/L Cya | 5-23 | 100 | 931223 |
| Fluoride* | 0.1–2.0 mg/L F ⁻ | 5-27 | 150 | 931227 |
| Iron 1* | 0.04–2.00 mg/L Fe | 5-25 | 150 | 931225 |
| Iron 2 | 0.04–2.00 mg/L Fe | 5-26 | 150 | 931226 |
| ■ pH 6.0–8.2 | 6.10-8.40 | 5-70 | 150 | 931270 |
| ■ Silica HR 200* | 10–200 mg/L SiO ₂ | 5-34 | 100 | 931234 |
| VISOCOLOR [®] Powder Pillows | | | | |
| Chlorine, free | 0.03–6.00 mg/L Cl ₂ | 7-20 | 100 1000 | 936220 936220.1 |
| Chlorine, total | 0.03–6.00 mg/L Cl ₂ | 7-21 | 100 1000 | 936221 936221.1 |
| оН | pH 6.2–8.2 | 7-22 | 100 | 936222 |
| Silica LR* | 0.02–2.10 mg/L SiO ₂ | 7-24 | 100 | 936224 |
| Silica HR* | 2–210 mg/L SiO ₂ | 7-25 | 100 | 936225 |
| Sulfate* | 15–200 mg/L SO ₄ ^{2–} | 7-23 | 100 | 936223 |
| NANOCOLOR [®] tests | | | | |
| Chlorine/Ozone | 0.05–2.50 mg/L Cl ₂ / 0.05–2.00 mg/L O ₃ | 0-17 | 20 | 985017 |
| ■ pH 6.5–8.2 | 6.10-8.40 | 0-72 | 100 | 91872 |



Ordering information

| Description | REF |
|---|--------|
| Compact photometer PF-3 Drinking Water, in a cardboard box, incl. manual, batteries and certificate | 919343 |
| Compact photometer PF-3 Drinking Water, in a case with foam inlay, incl. manual, batteries and certificate | 934402 |
| Reagent case with photometer PF-3 Drinking Water, | 934124 |
| with VISOCOLOR [®] ECO reagents* chlorine 2 (liquid), free and total, pH 6.0–8.2, fluoride, iron 2, chlorine dioxide | |
| Reagent case with photometer PF-3 Drinking Water, | 934125 |
| with VISOCOLOR® ECO reagents* chlorine 6 (solid), free and total, pH 6.0–8.2, fluoride, iron 2, chlorine dioxide | |

For soil analysis

Monitoring and reliable analysis is a corner stone in supporting and maintaining healthy, productive and biologically active soil. To effectively and efficiently plan all measures that affect the soil (fertilization, liming, etc.) it is crucial to determine the important soil parameters (nitrate, phosphate, potassium and ammonium). The PF-3 Soil in combination with the proven *VISOCOLOR®* ECO and *NANOCOLOR®* tests are ideally suited for mobile analysis of the soil quality right at the point of interest.



Tests suitable for the PF-3 Soil

| Tests | Measurement range | Test number | No. of tests | REF |
|----------------------------------|------------------------------------|-------------|--------------|--------|
| VISOCOLOR [®] ECO tests | | | | |
| Ammonium 3* | 0.2–4.0 mg/kg NH ₄ -N | 5-08 | 50 | 931208 |
| Potassium* | 40–300 mg/kg K ⁺ | 5-32 | 60 | 931232 |
| Nitrate | 2–28 mg/kg NO ₃ -N | 5-41 | 110 | 931241 |
| Phosphate* | 20–500 mg/kg PO ₄ -P | 5-84 | 80 | 931284 |
| NANOCOLOR® tests | | | | |
| Ammonium 3* | 0.08–4.60 mg/kg NH ₄ -N | 0-03 | 20 | 985003 |
| Ammonium 10* | 0.4–16.0 mg/kg NH ₄ -N | 0-04 | 20 | 985004 |
| Ammonium 50* | 2.0–80.0 mg/kg NH ₄ -N | 0-05 | 20 | 985005 |
| Potassium 50* | 40–1000 mg/kg K ⁺ | 0-45 | 20 | 985045 |
| ■ Nitrate 50* | 1–44 mg/kg NO ₃ -N | 0-64 | 20 | 985064 |
| Phosphate 5* | 4–100 mg/kg PO ₄ -P | 0-81 | 20 | 985081 |
| Phosphate 15* | 6–300 mg/kg PO ₄ -P | 0-80 | 20 | 985080 |



Ordering information

| Description | REF |
|--|--------|
| Compact photometer PF-3 Soil, in a cardboard box, incl. manual, batteries and certificate | 919341 |
| Compact photometer PF-3 Pool, in a case with foam inlay, incl. manual, batteries and certificate | 934202 |
| Reagent case with photometer PF-3 Soil, with VISOCOLOR [®] ECO reagents* nitrate, potassium, ammonium 3, phosphate | 934210 |
| Reagent case for soil analysis with photometer PF-3 Soil, incl. reagents and accessories for sample preparation* | 934220 |

For fishing waters

A regular water monitoring is the precondition for a healthy and productive living space for fishes. To avoid deficits and fish mortality, it is important to determine all relevant parameters and react promptly.

With it's three wavelengths, the compact photometer PF-3 Fish is perfectly suited for the analysis of fishing waters. Together with the proven *VISOCOLOR®* ECO and *NANOCOLOR®* tests, the PF-3 Fish is the perfect tool for a fast and reliable statement on the water quality.



Tests suitable for the PF-3 Fish

| Tests | Measurement range | Test number | No. of tests | REF |
|-------------------------------------|-----------------------------------|-------------|--------------|--------|
| VISOCOLOR [®] ECO tests | | | | |
| Ammonium 3* | 0.1–2.0 mg/L NH ₄ -N | 5-08 | 50 | 931208 |
| Chlorine 6 free, total and combined | 0.05–6.00 mg/L Cl ₂ | 5-17 | 200 | 931217 |
| Iron 1* | 0.04–2.00 mg/L Fe | 5-25 | 150 | 931225 |
| Iron 2 | 0.04–2.00 mg/L Fe | 5-26 | 150 | 931226 |
| ■ Silica* | 0.2–3.0 mg/L SiO ₂ | 5-33 | 80 | 931233 |
| Copper | 0.1–5.0 mg/L Cu ²⁺ | 5-37 | 100 | 931237 |
| Nitrate | 1.0–14.0 mg/L NO ₃ -N | 5-41 | 110 | 931241 |
| Nitrite | 0.01–0.15 mg/L NO ₂ -N | 5-44 | 120 | 931244 |
| ■ pH 6.0-8.2 | 6.10–8.40 pH | 5-70 | 150 | 931270 |
| Phosphate* | 0.2–5.0 mg/L PO ₄ -P | 5-84 | 80 | 931284 |
| Oxygen* | 1–8 mg/L O ₂ | 5-88 | 50 | 931288 |
| NANOCOLOR [®] tests | | | | |
| Ammonium 3* | 0.04–2.30 mg/L NH ₄ -N | 0-03 | 20 | 985003 |



Ordering information

| Description | REF |
|--|--------|
| Compact photometer PF-3 Fish, in a cardboard box, incl. manual, batteries and certificate | 919345 |
| Compact photometer PF-3 Fish, in a case with foam inlay, incl. manual, batteries and certificate | 934602 |
| Reagent case with compact photometer PF-3 Fish, incl. VISOCOLOR® ECO reagents, VISOCOLOR® HE reagents, QUANTOFIX® test strips and manual, batteries, certificate and accessories* | 934127 |

For COD analysis

COD levels play a major role in assessing communal and industrial waste water. As a sum parameter indicating the organic load of waste water, COD is one of the most important parameters in waste water analysis. Among others, COD values are used to calculate and monitor the purification capacities of waste water plants. The PF-3 COD together with the user friendly *NANOCOLOR*[®] tests is ideally suited to obtain information on the COD content (2–60000 mg/L) easily and directly at the place of sampling.



Tests suitable for the PF-3 COD

| Tests | | Measurement range | Test number | No. of tests | REF |
|-----------------|-----------|-------------------------------|-------------|--------------|--------|
| NANOCOLOR® test | | | | | |
| COD 40* | ISO 15705 | 2–40 mg/L O ₂ | 0–27 | 20 | 985027 |
| COD 60* | ISO 15705 | 5–60 mg/L O ₂ | 0–22 | 20 | 985022 |
| COD 160* | ISO 15705 | 15–160 mg/L O ₂ | 0–26 | 20 | 985026 |
| COD 600* | ISO 15705 | 50–600 mg/L O ₂ | 0–30 | 20 | 985030 |
| COD 1500* | ISO 15705 | 100–1500 mg/L O ₂ | 0–29 | 20 | 985029 |
| COD LR 150* | ISO 15705 | 3–150 mg/L O ₂ | 0–36 | 20 | 985036 |
| COD HR 1500* | ISO 15705 | 20–1500 mg/L O ₂ | 0-38 | 20 | 985038 |
| COD 4000* | | 400–4000 mg/L O ₂ | 0–11 | 20 | 985011 |
| COD 10000* | | 1.00–10.00 g/L O ₂ | 0–23 | 20 | 985023 |
| COD 15000* | | 1.0–15.0 g/L O ₂ | 0–28 | 20 | 985028 |
| COD 60000* | | 5.0–60.0 g/L O ₂ | 0–12 | 20 | 985012 |



Ordering information

| Description | REF |
|---|--------|
| Compact photometer PF-3 COD, in a box, incl. manual, batteries and certificate | 919342 |
| Compact photometer PF-3 COD, in a case with foam inlay, incl. manual, batteries and certificate | 934302 |
| NANOCOLOR® analysis case with compact photometer PF-3 COD, | 919212 |
| with additional space for 1 heating block NANOCOLOR $^{\oplus}$ VARIO C2 or NANOCOLOR $^{\oplus}$ VARIO Mini, 2 pipettes, | |
| 3 NANOCOLOR [®] tube tests and accessories | |

VISOCOLOR® ECO and NANOCOLOR® tests suitable for the PF-3

| VISOCOLOR® ECO tests | S | Measurement range | Test number | Version | No. of tests | REF |
|---------------------------------|-----------|--|-------------|---------|--------------|----------|
| Alkalinity TA | | 5–250 mg/L CaCO $_3$ | 5-04 | | 100 | 931204 |
| Ammonium 3* | | 0.1–2.5 mg/L NH4 ^{+ 1)} | 5-08 | | 50 | 931208 |
| Bromine | | 0.10–13.00 mg/L Br ₂ | 5-11 | | 200 | 931211 |
| Chlorine 1, free, total and | combined | 0.05–2.00 mg/L Cl ₂ | 5-35 | | 150 | 931235 |
| Chlorine 2, free, total and | combined* | 0.05–2.00 mg/L Cl ₂ | 5-15 | | 150 | 931215 |
| free Chlorine 2* | | 0.05–2.00 mg/L Cl ₂ | 5-16 | | | 931216 |
| Chlorine 6, free, total and | combined | 0.05–6.00 mg/L Cl ₂ | 5-17 | | 200 | 913217 |
| free Chlorine 6 | | 0.05–6.00 mg/L Cl ₂ | 5-19 | | 400 | 931219 |
| Chlorine dioxide* | | 0.10–3.80 mg/L CIO ₂ | 5-21 | | 150 | 931221 |
| Copper | | 0.1–5.0 mg/L Cu ²⁺ | 5-37 | | 100 | 931237 |
| Cyanuric acid* | | 10–100 mg/L Cya | 5-23 | | 100 | 931223 |
| Fluoride* | | 0.1–2.0 mg/L F⁻ | 5-27 | | 150 | 931227 |
| Iron 1* | | 0.04–2.00 mg/L Fe | 5-25 | | 150 | 931225 |
| Iron 2 | | 0.04–2.00 mg/L Fe | 5-26 | | 150 | 931226 |
| Nitrate | | 1.0–14.0 mg/L NO ₃ -N ¹⁾ | 5-41 | | 100 | 931241 |
| Nitrite | | 0.01–0.15 mg/L NO ₂ -N | 5-44 | | 120 | 931244 |
| Oxygen* | | 1–8 mg/L O ₂ | 5-88 | | 50 | 931288 |
| pH 6.0–8.2 | | 6.10-8.40 | 5-70 | | 150 | 931270 |
| Phosphate* | | 0.2–5.0 mg/L PO ₄ -P ¹⁾ | 5-84 | | 80 | 931284 |
| Potassium* | | 2–15 mg/L K ^{+ 1)} | 5-32 | • | 60 | 931232 |
| Silica* | | 0.2–3.0 mg/L SiO ₂ | 5-33 | | 80 | 931233 |
| Silica HR 200* | | 10–200 mg/L SiO ₂ | 5-34 | | 100 | 931234 |
| VISOCOLOR [®] Powder P | Pillows | | | | | |
| Chlorine, free | | 0.03–6.00 mg/L Cl ₂ | 7-20 | | 100 | 936220 |
| | | - | | | 1000 | 936220.1 |
| Chlorine, total | | 0.03–6.00 mg/L Cl ₂ | 7-21 | | 100 | 936221 |
| | | | | | 1000 | 936221.1 |
| pН | | pH 6.2–8.2 | 7-22 | | 100 | 936222 |
| Silica LR* | | 0.02–2.10 mg/L SiO ₂ | 7-24 | | 100 | 936224 |
| Silica HR* | | 2–210 mg/L SiO ₂ | 7-25 | | 100 | 936225 |
| Sulfate* | | 15–200 mg/L SO ₄ ^{2–} | 7-23 | | 100 | 936223 |
| NANOCOLOR® tests | | | | | | |
| Ammonium 3* | | 0.05–3.00 mg/L NH ₄ ^{+ 1)} | 0-03 | | 20 | 985003 |
| Ammonium 10* | | 0.20–10.0 mg/L NH ₄ ^{+ 1)} | 0-04 | • | 20 | 985004 |
| Ammonium 50* | | 1.0–50.0 mg/L NH ₄ ^{+ 1)} | 0-05 | • | 20 | 985005 |
| Chlorine/Ozone | | 0.05–2.50 mg/L Cl_2 / 0.05–2.00 mg/L O_3 | 0-17 | | 20 | 985017 |
| COD 40* | ISO 15705 | 2–40 mg/L O ₂ | 0-27 | | 20 | 985027 |
| COD 60* | ISO 15705 | 5–60 mg/L O ₂ | 0-22 | • | 20 | 985022 |
| COD 160* | ISO 15705 | 15–160 mg/L O ₂ | 0-26 | | 20 | 985026 |
| COD 600* | ISO 15705 | 50–600 mg/L O ₂ | 0-30 | | 20 | 985030 |
| COD 1500* | ISO 15705 | 100–1500 mg/L O ₂ | 0-29 | | 20 | 985029 |
| COD LR 150* | ISO 15705 | 3–150 mg/L O ₂ | 0–36 | | 20 | 985036 |
| COD HR 1500* | ISO 15705 | 20–1500 mg/L 02 | 0-38 | | 20 | 985038 |
| COD 4000* | | 400-4000 mg/L O ₂ | 0-11 | | 20 | 985011 |
| COD 10000* | | 1.00–10.00 g/L O ₂ | 0-23 | | 20 | 985023 |
| COD 15000* | | 1.0–15.0 g/L O ₂ | 0-28 | | 20 | 985028 |
| COD 60000* | | 5.0–60.0 g/L O ₂ | 0-12 | | 20 | 985012 |
| Nitrate 50* | | 0.3–22.0 mg/L NO ₃ -N ¹⁾ | 0-64 | | 20 | 985064 |
| pH 6.5–8.2 | | 6.10-8.40 | 0-72 | | 100 | 91872 |
| ortho and total Phosphate | ə 5* | 0.20–5.00 mg/L PO ₄ -P ¹⁾ | 0-81 | | 20 | 985081 |
| ortho and total Phosphate | | 0.30–15.00 mg/L PO ₄ -P ¹⁾ | 0-80 | | 20 | 985080 |
| Potassium 50* | | $2-50 \text{ mg/L K}^{+1)}$ | 0-45 | | 20 | 985045 |
| Polassium 50 | | | | | | |

* This product contains harmful substances which must be specially labeled as hazardous. For detailed information please see SDS.

¹⁾ For soil specific measurement ranges please refer to page 8 of the selection guide or www.mn-net.com/PF-3.

Additional tests and versions will be launched successively.

You can find all your current options at www.mn-net.com/PF-3.

Technical data

| PF-3 | |
|---------------------|---|
| Туре | LED photometer with microprocessor control, self-test and auto calibration |
| Optics | LED + interference filters |
| | Insensitive to external light for quick measurements without cuvette slot cover |
| Wavelengths | 3 wavelengths, depending on version |
| | Pool/Drinking Water: 450 nm/530 nm/590 nm Soil: 365 nm/450 nm/660 nm |
| | COD: 365 nm/450 nm/595 nm |
| | Fish: 450 nm / 530 nm / 660 nm |
| Wavelength accuracy | \pm 2 nm, bandwidth at half transmission 10 nm–12 nm |
| Light source | LED |
| Detector | Silicon photodiode |
| Measuring modes | NANOCOLOR [®] tube tests, VISOCOLOR [®] ECO tests and |
| | VISOCOLOR [®] Powder Pillows |
| Cuvette slot | Tubes 16 mm OD |
| Memory | 50 results |
| Display | Backlit graphic display, 128 x 64 pix, all important data at a glance: result with unit, date, time |
| Operation | Self-explanatory menu guidance, plastic foil keypad, test selection via parameter lists |
| Interface | Mini-USB |
| Update | Free of charge via internet / PC |
| Operating range | 10 °C–40 °C, up to 80 % relative humidity (non-condensing) |
| Power supply | 3 AA batteries, rechargeable batteries, USB interface; optional internal accu-pack |
| Housing | Shock-resistant, waterproof and dustproof, according to IP68 |
| Dimensions | 170 mm x 95 mm x 68 mm |
| Weight | 0.5 kg |
| Warranty | 2 years |
| CE | CE-certified |
| | |

www.mn-net.com

MN EN ISO 9001 CERTIFIED

MACHEREY-NAGEL GmbH & Co. KG \cdot Neumann-Neander-Str. 6–8 \cdot 52355 Düren \cdot Germany DE/International: Tel.: +49 24 21 969-0 Fax: +49 24 21 969-199 E-mail: info@mn-net.com

MACHEREY-NAGEL

CH: Tel.: +41 62 388 55 00 Fax: +41 62 388 55 05

FR: Tel.: +33 388 68 22 68 Fax: +33 388 51 76 88 E-mail: sales-ch@mn-net.com E-mail: sales-fr@mn-net.com

