
sartorius stedim
biotech

Directions for Use | Gebrauchsanleitung | Mode d'emploi
Instrucciones para el uso | Istruzioni per l'uso

Midisart[®] 2000

17804 | 17805 | 17877

82521-060-50

Directions for Use of Midisart 2000

17804 | 17805 | 17877

1. Product Description and Applications

Midisart 2000 ready-to-connect filtration units consist of a hydrophobic (PTFE) membrane filter in a polypropylene housing. They can be ordered sterile and non-sterile. The filtration unit can be repeatedly autoclaved. They are available with 0.2 µm and 0.45 µm pore size membranes, and are ideally suitable for all small-scale air|gas sterilizing purposes, such as:

- Filling and transfer vessels
- Holding and storage tanks for sterile media
- Autoclaves
- Small (single-use) fermenters
- Disposable bags | bottles and tubing systems
- Carboys
- Solvent filtration
- Single-use assemblies

The excellent chemical stability of these practical units also makes them useful for removal of particles from chemicals and organic solvents.

The female luer connection in the hose nipple allows the unit to be firmly fitted onto an appropriate syringe.

2. Order Numbers

Pore size in µm	Inlet and outlet connectors	Box of 12, sterile, individually packaged units	Box of 25, sterile, individually packaged units	Box of 100, bulk-packaged, non-sterile units	Plastic bucket of 100, bulk packaged, non-sterile units
0.2	Hose nipple	17805-E	17805-G	17805-UPN	17805-UQN
0.2	Modified hose nipple	-	-	17877-UPN	
0.45	Hose nipple	17804-E	17804-G		
0.2	1/8" NPT thread	17805-NPE	17805-NPG		
0.45	1/8" NPT thread	17804-NPE	17804-NPG		

Midisart units can also be used for the filtration of aqueous liquids, but the hydrophobic PTFE membrane filter must first be wetted with a water-miscible solvent, such as ethanol.

3. Specifications

Filter material: PTFE reinforced with polypropylene gauze

Housing material: Polypropylene

Filtration area: 20 cm²

Housing diameter: 64 mm

Priming volume: approx. 3 mL

Hold-up volume: approx. 1 mL before, 0.5 mL after bubble point

Max. operating pressure: 300 kPa (3 bar)

Water penetration point: approx. 300 kPa (0.45 µm), 400 kPa (0.2 µm)

Max. temperature (autoclave): 134°C

Presterilized with: Ethylene oxide

Biosafety: Materials pass USPXXII Plastics class VI

Bubble point (isopropanol, 60 %):

0.45 µm: ≥0.9 bar

0.2 µm: ≥1.1 bar

Connections: Inlet and outlet: Conical hose nipple for hoses with an inner diameter of 6–12 mm, modified hose nipple for hoses with an inner diameter of 4–12 mm, or

Inlet and outlet: 1/8" NPT thread

Typical Air Flow Rates per Pore Size

Differential Pressure	0.2 µm	0.45 µm
Δp = 0.02 bar	1.1 L/min	1.8 L/min
Δp = 0.05 bar	2.9 L/min	4.6 L/min
Δp = 0.10 bar	5.0 L/min	8.5 L/min

1 bar = 100 kPa = 14.5 psi

4. Chemical Compatibility

The materials (PP and PTFE) used in Midisart give it excellent compatibility with the products listed below. However, its compatibility can be affected by various factors, such as temperature, concentration, composition, etc. We therefore recommend that you test whether Midisart is compatible with the particular medium you wish to filter.

Media

Acetic acid (concentrated)

Acetone

Acetonitrile

n-Butanol

Cellosolve (ethyl)

Chloroform

Diethylacetamide

Dimethyl formamide

Dimethylsulfoxide

Dioxane

Ethyl acetate

Ethanol

Ethylene glycol

Freon TF

Gasoline

1 N hydrochloric acid

Hexane

Isobutanol

Isopropanol

Methanol

Methylene chloride

Methyl ethyl ether

Methyl ethyl ketone

Pentane

Sodium hydroxide (5%)

Tetrahydrofuran

Toluene

Trichloroacetic acid

Trichloroethane

Water

Xylene

In the interest of the further development of Sartorius Stedim Biotech products, we reserve the right to make changes without notice.

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Phone +49.551.308.0
Fax +49.551.308.32 89
www.sartorius-stedim.com

Copyright by
Sartorius Stedim Biotech GmbH,
Goettingen, Germany.
All rights reserved. No part
of this publication may
be reprinted or translated in
any form or by any means
without the prior written
permission of Sartorius Stedim
Biotech GmbH.

The status of the information,
specifications and illustrations
in this manual is indicated
by the date given below.
Sartorius Stedim Biotech GmbH
reserves the right to make
changes to the technology,
features, specifications and
design of the equipment
without notice.

Status:
April 2015,
Sartorius Stedim Biotech GmbH,
Goettingen, Germany

Specifications subject to change without
notice. Printed in the EU on paper
bleached without chlorine.
W1_Midisart 2000 · KT
Publication No.: SLD6017-p150402
Ver. 04 | 2015

List of Sartorius material numbers applying to EPA-FIFRA

17804	17805---NPG	J1-17805
17804-----E	17805---PCM	
17804-----G	17805---SUN	
17804-----M	17805---TCN	
17804---CPW	17805---UFQ	
17804---NAN	17805---UKN	
17804---NP	17805---UPB	
17804---NPE	17805---UPN	
17804---NPG	17805---UPQ	
17804---NQN	17805---UQ1	
17804---SUN	17805---UQN	
17805	17805---UZQ	
17805-----B	17806---ATJ	
17805-----E	17809---SAN	
17805-----G	17809---SUN	
17805-----M	17809---UNN	
17805---BME	1780D---UPN	17811---UNN
17805---CDE	17811---UNN	
17805---CPW	17812---UNN	
17805---CXE	17839---UNN	
17805---DCE	17877---AI1	
17805---NAN	17877---DRN	
17805---NP	17877---UPN	
17805---NPE	18072	