



CYAN CL002 Rotator

Suitable for mixing larger volumes or for agglutination tests.



Choose for **CYAN CL002** and get the following advantages:

- Economical and practical instrument
- Runs at steady and stable speed
- 100 RPM speed is indicated (advised speed for RPR Carbon tests)
- Compact
- Quiet
- Vibration free
- Includes variable speed and timer control
- Platform is furnished with anti slip pad
- Fixation spring for larger objects provided



Technical data:

- Speed for 0-1 kg: 20 - 240 RPM
- Speed for 1-2 kg: 20 - 150 RPM
- 0 - 30 minutes timer
- Platform dimensions: +/- 23 x 23 cm
- Instrument dimensions: +/- 26 x 29,5 x 14 cm (WxDxH)
- Load capacity: +/- 2 kg
- Orbit diameter: +/- 20 mm
- Gross weight: +/- 3,7 kg
- Power supply: 110 V - 60 Hz (CL002-110)
or 220 V - 50 Hz (CL002-220)

Order code:

- CL002-110: Rotator (110 V - 60 Hz)
- CL002-220: Rotator (220 V - 50 Hz)
- CL002-220/60: Rotator (220 V - 60 Hz)



Specification subject to change without prior notice.

Cypress Diagnostics is a quality oriented and ISO certified company with CE approved products
Address : Langdorpsesteenweg 160 3201 Langdorp Belgium

